

# **INTRODUCTION TO CUTEGEL**

# What is Cutegel?




**1** Non-animal Hyaluronic Acid (HA) filler

**2** HA - natural substance in human body

**3** Low probability of allergic reactions

**4** Class IV medical device

# Cutegel Features



# Strict Quality Control

## Safe Ingredient




EDQM certified  
HA ingredient

## Thorough Wash Out

No BDDE residue  
due to thorough  
wash out process


## Minimized Side Effect

Low endotoxin and  
protein residues



[GC Analysis: BDDE experiment result ]


# Certificates




Manufacturer



GMP Certificate



Venture Business



R&D Center



ISO14001



ISO9001



Cutegel (HA Dermal Filler)



Hibarry (Anti-adhesion Barrier Gel)



ISO13485



Patent

CE Marking

# Clinically Proven Efficacy & Safety

: Non-inferiority of Cutegel in Comparison to Restylane product based on WSRS and GAIS scores




Fig. \*WSRS scores after injection

\*\*GAIS scores comparison evaluated by investigator


|           | +8 weeks   | +16 weeks  | +24 weeks  |
|-----------|------------|------------|------------|
| Cutegel   | 1.6 ±0.63  | 1.59 ±0.74 | 1.37 ±0.57 |
| Restylane | 1.59 ±0.63 | 1.59 ±0.74 | 1.37 ±0.57 |

\*\*GAIS scores comparison evaluated by subjects


|           | +8 weeks   | +16 weeks  | +24 weeks  |
|-----------|------------|------------|------------|
| Cutegel   | 1.47 ±0.84 | 1.41 ±0.92 | 1.31 ±0.92 |
| Restylane | 1.53 ±0.78 | 1.42 ±0.90 | 1.19 ±0.98 |

# Even Particle Structure


Compared to other companies, Cutegel is consisted of even HA particles through optimized particle purification method.



CUTEGEL




Filler A




Filler B

# Graph of Particle Size



# Various Particle Size for Different Uses

Cutegel provides different products with various HA particles that are appropriate for different layers.



Cutegel Aqua



Cutegel S



Cutegel N



Cutegel Max



Cutegel Max 1400

# CUTEGEL Line-up

the ageless fantasy C U T E G E L

## CUTEGEL Aqua S

0.9ml  
20mg/ml  
Upper Dermis



## CUTEGEL Aqua

2.2ml  
20mg/ml  
Upper Dermis



## CUTEGEL S

0.7 ml  
20mg/ml  
Upper to Mid Dermis



## CUTEGEL N

1.0 ml  
20mg/ml  
Mid to Deep Dermis



## CUTEGEL Max

1.1ml  
20mg/ml  
Subcutaneous Tissue



## CUTEGEL Max 1400

1.1ml x 2 Syringe  
20mg/ml  
Subcutaneous Tissue



## CUTEGEL Ultra

2.0ml  
20mg/ml  
Subcutaneous Tissue



## CUTEGEL Strong

3.0ml  
20mg/ml  
Subcutaneous Tissue




# CUTEGEL Indications



CUTEGEL N



CUTEGEL MAX



CUTEGEL N



CUTEGEL N



CUTEGEL MAX



CUTEGEL AQUA



CUTEGEL MAX

CUTE-PATCH



CUTEGEL N



CUTEGEL S

**Skin  
Hydration**



CUTEGEL AQUA

# Long-lasting Hydrolifting Injection

## CUTEGEL AQUA



- Longer duration due to cross-linked HA
- For whole face, décolté , cleavage, and etc
- Effective for moisturization and pore tightening
- Firmer skin with better texture, less surface roughness, and more elastic
- 3 times every 8 weeks recommended
- 1 injection effective for more than 6 months

# Be More Attractive With CUTEGEL S



- Each particle size 200um
- For around mouth and fine wrinkles
- 30G needle enables smooth injection
- Specialized for lip lines
- Lasts for 6 to 12 months

# Multiplayer Filler CUTEGEL N



- Each particle size 600um
- 27G needle fit for midsize particles
- Used most extensively in various areas
- Mainly for nasolabial fold and nose
- Naturally restores face volume
- Lasts for 6 to 12 months

# Natural Volume Filler CUTEGEL MAX



- Each particle size 1000um, 1400um
- Volume filler with largest particle size
- Specialized for volumizing
- For forehead, cheek, chin, and etc
- Natural volumizing resulting in increased satisfaction
- Lasts for longer than 12 months

# Natural Volume Filler CUTEGEL Ultra & Strong



**2 ml**



**3 ml**

- Each particle size 1400um
- Volume filler with largest particle size
- Specialized for Body
- Natural volumizing resulting in increased satisfaction
- Lasts for longer than 12 months

**THANK  
YOU**